Introduction

We would like to thank you for your decision to purchase this HAAG-STREIT product.

If the instructions in this manual are carefully followed we are confident that this product will give you reliable and trouble-free usage.

We will be forced to disclaim any warranty and liability if the instrument is altered in any way or if routine maintenance is neglected or not carried out according to factory specifications.

Purpose of use

The **contact glasses** serve in diagnosis and / or laser treatment of the human eye.

They are usually used, at room temperature, in consulting rooms, clinics, hospitals or teaching facilities.

They are used by Ophthalmologists and also Optometrists, if they are qualified by their respective authorities to do so.

NOTICE

All rights to technical alterations reserved.

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	HERA INTERNATIONAL	
1	Contents	page
2	Operating safety	5
3	Description of the product	11

Contact glasses for diagnostics		
3.1	901	Contact glass 13
3.2	902	1 mirror contact glass 15
	902 S	1 mirror contact glass 15
3.3	905	2 mirror contact glass 17
	905 S	2 mirror contact glass 17
3.4	904	1 mirror contact glass 19
3.5	General information 3 mirror gl. 21	
3.6	903	3 mirror contact glass 21
	903 S	3 mirror contact glass 21
3.7	630	3 mirror contact glass 23
3.8	906	3 mirror contact glass 23
	906 S	3 mirror contact glass 23
	907	3 mirror contact glass 23
3.9	1210	Endothelium contact glass 25
Contact glasses for Laser application		
3.10	CGAL	Gonioscopic contact gl. 29
3.11	903 L	3 mirror contact glass 31
	630 L	3 mirror contact glass 31
	906 L	3 mirror contact glass 31
	907 L	3 mirror contact glass 31
3.12	CGRL	Retina contact glass 33
3.13	CGIL	Irdectomy contact glass 35
3.14	CGVL	Vitrectomy contact glass 37
3.15	CGPL	Capulotomy contact gl. 39
3.16	RETINA 145 L	Panfundus cont. gl. 41
4	Operating the product	
A	Maintenance	
A.1	Cleaning and disinfection 45	
A.2	Accessories disinfection / storing 45	
Short How to disinfect? 47		
B	Technical specifications	
		49

Scleraler Rand	Bord scléral	Scleral flange	
Scleraler Rand	Bord scléral	Scleral flange	
Scleraler Rand	Bord scléral	Scleral flange	
Kurze Distanz	Courte distance	Short distance	
Neugeborene	Nouveau-nés	Newborn	
+ Scleraler Rand	+ Bord scléral	+ Scleral flange	
Kinder -4 Jahre	Enfants -4 ans	Children -4 yea	
Kurze Distanz	Courte distance	Short distance	
Neugeborene	Nouveau-nés	Newborn	
Kinder -4 Jahre	Enfants -4 ans	Children -4 yea	

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2 Operating safety

2.1 Safety Regulations

The following instructions should be observed to ensure the safety of patients:

Avoid examinations in case of eye infections or injured corneas!

- There is no contraindication for applications of contact glasses. Medical assessment is necessary and due care is to be exercised in case of infections of the eye and injuries to the cornea.

Only clean and disinfected contact glasses should be used!

- Examinations should only be conducted with cleaned and disinfected contact glasses.
- Contact glasses are not dispatched in a disinfected state and should always be disinfected before use.
- The disinfecting of contact glasses eg. with cotton wool or damp cloth without putting them into disinfecting fluid is insufficient and not allowed.
- Incorrect medical disinfecting can cause patients to be infected and can also cause damage to the eyes.
- Residue of the disinfecting may cause provoking and caustic reactions on the patient's eye, therefore rinse thoroughly.
- Cleaning and disinfecting is to be carried out according to Appendix A and short manual.

No damaged contact glasses should be used!

- Defective contact glasses should not be used, in order to avoid any damage to the cornea of the patient.

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Before each use the contact surface of all contact glasses should be checked for contamination or damage (such as scratching, rips or sharp edges). The slit lamp microscope is best used for this purpose at a magnification of 10 x to 16 x.

• The presence of cracks or defective filling on the front of the contact glass can allow disinfectant to seep into the hollow interior, leading to eye irritation for the patient. Contact glasses that have allowed condensation to enter may no longer be used.

2.2 Transportation, storage and installation

• The contact glasses may only be exposed to the following environmental/ambient conditions for transportation- and storage purposes (opened or in the packaging) for a maximum of 15 weeks:

a range of ambient temperatures
from -10 °C to +60 °C

a range of air pressures
from 500 hPa to 1060 hPa

a relative humidity range from 10 % to 90 %, whereby the contact glasses should be allowed to adjust to the ambient temperature of the new location for several hours in the packaging (to avoid the forming of condensation).

• The contact glasses should be checked for damage after unpacking.

• Defective equipment should always be returned in an appropriate packaging.

2.3 Notes on usage

• Only qualified and trained personnel may operate and clean the contact glasses.

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- The training of the operating personnel is the responsibility of the product user.
- In order to avoid possible eye damage, the sun should never be viewed through the ocular.
- The contact glasses must be checked after exposure to external force (e.g. unintentional shocks or being allowed to fall) and if necessary they are to be returned to the manufacturer for repair.
- The manufacturer will not be liable for any loss or damage resulting from unauthorised intervention. All ensuing warranty claims will be deemed null and void as a result of unauthorised interaction.

2.4 Statutory requirements

- The contact glasses have been designed and constructed to conform with **EN ISO 13485** and **EN 30993-3** standards. Manufacturing procedures, testing, commissioning, maintenance and repair are conducted under the observance of Swiss and International regulations.

- The contact glasses are designated as Class I under the EEC Directive 93/42 for medical equipment products. The CE designation was then granted on the basis of this Directive as per Conformity Module A.

- Statutory accident regulations are also to be observed.

2.5 Caution

- The safety regulations displayed in the operating instructions are to be observed with special care.

CAUTION Strictly observe all warning notices!

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3 Description of the product

The contact glasses, produced in conjunction with Prof. Goldmann, for the examination of the angle of the iris and the ocular fundus of the eye with the slit lamp are valuable diagnostic accessories, which have proved their use in ophthalmic practice and research.

In conjunction with the slit lamps of the 900® range, they give a binocular and stereoscopic view of the optical section of the entire eye, even under unfavorable conditions (cloudy media, small pupils, myopia, etc.).

The handling of the contact glasses is enhanced, whilst working with the slit lamp 900®, by the good working distance between the patient and the examiner.

Designation

- The contact glasses designated with an **"S"** possess a special scleral flange, which helps to create a better interface between the contact glass and the eyelid.
- Laser glasses are designated with an **"L"**.

Material


- Contact glasses only for **diagnostic purposes** are made of acrylic glass.
- Contact glasses for **laser treatment** applications consist of mineral glass, or have a coated mineral protection glass. The coating for antireflection is reached with a highly loadable coating and with high adhesive strength for the applicable laser wavelengths.

Stery Cup

- A symbol indicates which contact glasses can be used with HAAG-STREIT Stery Cup.

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Direct glasses for diagnostic purposes

3.1 Contact glass 901

For the binocular slit observation of the pupil, the macula, its surroundings out to 30° as well as the central vitreous body sections.

The refractive power is - 64 D. The linear, axial and angular enlargement does not depend on the refraction of the eye to be examined, which is of importance for measurements taken at the ocular fundus and in the vitreous body.


The contact glass has a weight of 1.5 grams and consists of two parts: the optically effective body of the glass with the haptic and the funnel to facilitate the handling of the contact glass. The effective optical surface has a diameter of 12 mm. Optic made of acrylic glass.

Examination of the ocular fundus


Before an examination with the fundus contact glass, a maximum mydriasis is desirable. The contact glass is positioned under topical anesthesia, using Methocel 2 % (see under 4 "Operating the product").

The patient's eye is brought into the required position by means of the fixation device and the area to be examined observed at 10 to 16 x magnification using the brightest possible, narrow slit.

During examination, the binocular, stereoscopic slit observation procedure is primarily to be followed, for which purpose the widest possible angle between the illumination and the microscope is of advantage. At the commencement of the examination, the beam is small and then enlarged, as soon as the slit image can be viewed binocularly.

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Mirror contact glasses for diagnostic purposes

3.2 One mirror contact glass 902 / 902 S

For the examination of the anatomical conditions in the area of the front of the iris.

This lens weighs 4.2 grams. Optic made of acrylic glass.

The expanse and form of the iridocorneal angle can only be ascertained by employing a narrow slit beam.

The slit lamp 900®, with its swivelling and declinable illumination, is primarily capable of providing an assessment of the expanse of the entire angle of the iris, which is of special importance in glaucoma diagnostics.


Examination of the angle of the iris

A magnification of 10 x is best for initial adjustments and 16 x during the actual examination.

Illumination of the iridocorneal angle is effected with a narrow slit via the contact glass mirror. An angle of around 10° is set between the illumination and the microscope (with the slit lamp 900® in the lateral dwell position). When an angle greater than 15° is employed between the microscope and the illumination, then the slit image will not be correctly focussed at the normal setting of the ocular.

An unfocussed slit at a degree less than 15° can be corrected by altering the position of the iridocorneal angle lens.


Examination of the lateral iridocorneal sections under diffused illumination conditions of the sclera sinus is also possible. For this, the illumination (after slackening the centering screw on the slit lamp 900®) is to be moved away from the central position.

• Gonio
Einspiegel-Kontaktglas 902
Verre de contact à un miroir 902
One mirror contact glass 902

• Gonio
Einspiegel-Kontaktglas 902 S
Verre de contact à un miroir 902 S
One mirror contact glass 902 S

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Examen de l'angle iridocornéen

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3.3 Two mirror contact glass 905 / 905 S


Two mirrored surfaces, which each incline at an angle of 62° towards the front surface, have been cut into the acrylic glass cone. This makes it possible to examine the opposite iridocorneal angle without turning the contact glass.

The expanse and form of the angle of the iris can be correctly ascertained by using a narrow slit on the slit lamp.

The slit lamp 900®, with its swivelling and declinable illumination, is primarily capable of providing an assessment of the expanse of the entire angle of the iris, which is of special importance in glaucoma diagnostics.

Examination of the angle of the iris

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3.4 One mirror contact glass 904


This one mirror contact glass, made of acrylic glass, serves to examine the ora serrata, the pars plana and the extreme peripheral parts of the vitreous body and the ocular fundus.

Examination of the ora serrata


When the eye has been anesthetized, 2 drops of Methocel 2 % are to be inserted into the concave part of the contact glass. It is of advantage to press down the lower lid by means of the depressor, have the patient look upwards and then draw the upper lid upwards with the thumb of the other hand. The ora glass is then placed on the cornea and, at the same time, the patient requested to look straight to the front.

A dilation of the pupil by more than 8 mm will enable observation of the ora serrata and its vicinity.

For cleansing and disinfection purposes, the ora glass is to be dismantled into its three component parts.

• Ora serrata
Einspiegel-Kontaktglas 904
Verre de contact à un miroir 904
One mirror contact glass 904

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3.5 General information on three mirror contact glasses

Examination of the entire ocular fundus and the iridocorneal angle

The advantage of a longer mirror is that it often permits binocular observation of the lateral sections of the ocular fundus.

Zone 1 = Lens (1)
Examination of the ocular fundus in the 30° zone, mostly with a magnification of 10 x

Zone 2 = Mirror (2) with a declination of 73°
Observation of the zone outside of the 30° range


Zone 3 = Mirror (3) with a declination of 66°
Observation of the peripheral sections of the ocular fundus, and under favorable conditions, also of the ora serrata

Zone 4 = Mirror (4) with a declination of 59°
Observation of the vitreous body and ocular fundus sections neighboring the ora serrata as well as a gonioscopic examination

3.6 Three mirror contact glass 903 / 903 S

The classic contact glass from Professor Goldmann!

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

Mirror positioning 1 - 4 and zones in the eye

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

3.7 Three mirror contact glass 630

This contact glass is particularly suitable for use with slit lamps which have a smaller distance between the eye and the reduction prism than the HAAG-STREIT 900®.

• Fundus tridoccornea
(Baby)
Dreispiegel-Kontaktglas 630
Verre de contact à 3 miroirs 630
Three mirror contact glass 630

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
3.8 Three mirror contact glasses for children 906 / 906 S and 907


Newborn: 906 / 906 S
Up to 4 years old: 907

The mirrors of these three mirror contact glasses, especially designed for children, are smaller than the normal three mirror contact glasses, but the declinations of the mirrors are identical.



Two models are available. One with an internal diameter of the cornea part of 10 mm, for newly born children (model 906) and the other with an 11 mm diameter for children up to 4 years of age (model 907).

• Fundus tridoccornea
(Baby)
Dreispiegel-Kontaktglas 906
Verre de contact à 3 miroirs 906
Three mirror contact glass 906

• Fundus tridoccornea
(Baby)
Dreispiegel-Kontaktglas 906 S
Verre de contact à 3 miroirs 906 S
Three mirror contact glass 906 S

• Fundus tridoccornea
(Children)
Dreispiegel-Kontaktglas 907
Verre de contact à 3 miroirs 907
Three mirror contact glass 907

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3.9 Contact glass 1210 endothelium according to Eisner

For the examination of the corneal endothelium and epithelium, magnification 2,2 x

Using the contact glass

- After removing the protective plug (4) fill the liquid reservoir with Methocel 2 %
- Insert into lid gap with release opening (1) resp. **marking (2) facing upwards**

Locating specular reflections

- Adjust slit lamp
 - Magnification:
 - Slit width:
 - Angle between slit lamp and microscope:
- Focus optical system with silver-coloured ring (3) and centre up
- Slide the slit lamp until gradually appear in sharp focus: the reflection of the optics (5), the corneal areas to be examined, the epithelium (6) respectively the endothelium (7)
- Adjust the contact glasses until mirror reflections appear illuminated

Examine specular reflections



- Slit lamp adjustment
- Magnification:
 - Slit width:
 - Angle between slit lamp and microscope:
- Examine periphery of cornea
- Guide patient's line of vision in desired direction
 - Adjust contact glass

Disinfecting

Insert protective plug (4) into contact glass. Now proceed as with other contact glasses.

• Endothelium
Kontaktglas nach Eisner 1210
Verre de contact d'après Eisner 1210
Contact glass according to Eisner 1210

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Use of the eyepiece with comparison grid according to McIntyre

Eyepiece 25 x for slit lamp BM 900®

- Insert eyepiece half-way into microscope and hold against cylinder housing with one hand.
- Use your other hand to loosen the knurled ocular refraction ring (2) by turning it anticlockwise through about 10°.
- De-focus the comparison grid by turning the smooth occluder (1) anticlockwise. The knurled ocular refraction ring should not turn at this point.
- Slowly turn the occluder clockwise, until the comparison grid is again sharply focussed.
- Secure the knurled ocular refraction ring.
- Push the eyepiece fully into the barrel.
- Align centre-marks vertically.
- Use knurled ocular refraction ring on **microscope** to adjust dioptre setting to zero.

Eyepiece 12.5 x for slit lamp BQ 900®

When adjusting the eyepiece to the user's refraction, look towards a light surface and adjust the eyepiece from the + side by turning the knurled ocular refraction ring with dioptre scale (3) until the eyepiece crosshairs appear sharply focussed.

Estimating the cell density of the endothelium

without contact glass

- Microscope magnification 40 x.
- Reading taken in column 40 x.

with contact glass 1210 (Eisner)

- Microscope magnification 25 x (overall magnification: microscope magnification x contact glass magnification = 55 x).
- Take reading of control sample in column 2.2 x 25 x

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Mirror contact glasses for laser treatment applications

3.10 Gonioscopic contact glass CGAL

The CGAL gonioscopic contact glass was designed by Roussel and Fankhauser for laser treatments of the chamber angle. It is made of laser resistant glass and is antireflection coated¹⁾. It reduces the laser spot by a factor of 1.5 as compared to a Goldmann lens. Therefore, lower energies are needed for cutting strands in the cataract incision with the YAG laser. Thanks to its magnifying effect, the CGAL contact glass also improves aiming accuracy and success rate in laser trabeculoplasty.

In office examination, the CGAL contact glass allows the visualization of the fine structures of the trabeculum with high resolution and produces brilliant photographs and video recording due to the low reflectivity coating in the visible spectral region.

Optical Properties

The convex entrance surface preserves the laser beam only convergence without introducing aberrations. It is this combination that results in minimal laser spot size and high optical resolution.

Magnification 1.5 x
Mirror angle 58°

¹⁾ Antireflection coating for visual observation and laser treatment:
450 - 650 nm and 1064 nm

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
The aiming point is the center of curvature C of the spectral entrance surface and as such, free of spherical aberrations and coma.

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
The beam aberrations (θ) are significantly lower with the CGAL than with any Goldmann types (γ: tilting angle of the contact glasses).

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



3.11 Three mirror contact glasses
for laser treatment applications

See description in 3.5 to 3.8

Acrylic lens with lumenized-coated mi-
neral protective glass for laser treatment¹⁾.

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Direct contact glasses for laser treatment applications

3.12 Retina contact glass CGRL

The CGRL is a wide angle contact glass designed for laser photocoagulation and diagnosis of the retina. It is made of laser resistant glass and is antireflection coated¹⁾. The optical design with a concave entrance surface allows a wide angle observation and laser treatment without intermediate imaging and without the associated intermediate focal spot inherent to indirect ophthalmoscopy.

The result is a small, lightweight contact glass which increases the field of view, without increasing the working distance, while providing an excellent image quality and good resolution. The whole posterior pole including the temporal vascular arcades, can be viewed simultaneously. It is used for axial and paraxial photocoagulation.

Optical Properties

The CGRL Retina contact glass creates an upright virtual image of the fundus. The cone angles of both laser and observation beam are reduced from 6° in air to 3,2° in vitreous. As a consequence the field of view is increased from 46° to 64°, which allows observation of the whole posterior pole.


Dioptric power in air	- 50 D
Magnification ²⁾	0,75 x
Field of view	64°
Spot diameter on the (Retina for 50 µm in air)	70 µm

¹⁾ Antireflection coating for visual observation and laser treatment: 450 - 650 nm


²⁾ The magnification with the Goldmann lens is arbitrarily set equal to 1.

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3.13 Iridectomy contact glass CGIL

The CGIL iridectomy contact glass was designed by Riquin, Fankhauser et al. to perform full thickness peripheral laser iridectomies. The use of this contact glass increases safety and efficacy of the laser procedure. It increases the energy density on the iris while minimizing it on anterior and posterior structures, such as cornea and lens capsule. The magnifying effect allows for improved aiming accuracy and higher success rate. Perfect handling is facilitated by the large diameter. It is made of laser resistant glass and is antireflection coated¹⁾.

Optical Properties

The CGIL contact glass increases the cone angle of the laser beam and therefore it increases the beam diameter on out-of-focus structures, while decreases the focal spot. The effective focal spot reduction is achieved by minimizing spherical aberrations and coma through the application of the Young-Weierstrass theorem: The lens has been designed with the applanatic point of the entrance surface located on the iris, at three millimeters posterior to the cornea.

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3.14 Vitrectomy Contact Glass CGVL

The CGVL vitrectomy contact glass was designed by Rol, Fankhauser et al. for photodisruptive YAG laser procedures in the posterior vitreous. While structures in the anterior to mid-vitreous may be treated with the CGPL or without any contact glass, safety and efficacy of photodisruption are increased in the deeper vitreous with the CGVL contact glass. Its magnifying effect and the possibility to visualize the retina allow for improved aiming accuracy. It is made of laser resistant glass and is antireflection coated¹⁾.

Optical Properties

The CGVL contact glass increases the cone angle of the laser beam and therefore it increases the beam diameter on out-of-focus structures, while it decreases the focal spot. The effective focal spot reduction is achieved by minimizing spherical aberrations and coma: The center of curvature of the entrance surface - a point free of spherical aberrations and coma - coincides with the center of curvature of the retina, while the other aplanatic point of the entrance surface is located on the retina.

Dioptric power in air - 33 D
Cone angle magnification 1 - 1.4 x

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3.15 Capsulotomy Contact Glass CGPL

The CGPL capsulotomy contact glass was designed by Riquin, Fankhauser et al. for the dissection of opacified posterior lens capsules and membranes in the pupillary and retropupillary space with the YAG laser. This contact glass increases safety and efficacy of the laser procedures. It lowers the minimal laser energy that is necessary for disruption of the capsule and it reduces the likelihood of IOL pitting. The magnifying effect allows for improved aiming accuracy, which is particularly important in the presence of a lens implant. Beam handling is facilitated by the large diameter of the contact glass. It is made of laser resistant glass and is antireflection coated¹⁾.

Optical Properties


The CGPL contact glass increases the cone angle of the laser beam and therefore it increases the beam diameter on out-of-focus structures, while it decreases the focal spot. The effective focal spot reduction is achieved by minimizing spherical aberrations and coma through the application of the Young-Weierstrass theorem: The contact glass has been designed with the applanatic point of the entrance surface located in the retropupillary space, at eight millimeters posterior to the cornea.

Dioptric power in air + 9 D
Cone angle magnification 1.5 x

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
Optical configuration of the CGPL capsulotomy contact glass with a pseudophakic eye.

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3.16 Contact glass RETINA 145 L

The Retina 145 L direct contact glass is a panfundus wide-angle viewing device designed to facilitate the problem-free diagnosis and treatment of the retina as far as the equator.

With an image magnification of 0,7 x, this direct contact glass provides highly detailed views of even subtle changes in the fundus.

The quality of the laser-scan image is as high on the periphery as it is in the centre. There is no major image-distortion.

Simultaneous visualization of the posterior pole and centre periphery permits optimum orientation at any time at the back of the eye, thus minimizing the risk of undesired laser coagulations in the area of the posterior pole.


In short, the Retina 145 L panfundus direct contact glass is an ideal tool for unlimited use in diagnosis and laser therapy in the area of the retina, including the equator. Its suitability for use in the diagnostic and therapeutic area makes this direct contact glass available for the entire range of retinological applications.

Dioptric performance in air	86.5 D
Magnification	- 0.7 x
Field of view	145°
Laser-point diameter	70 µm
(on the retina for 50 µm in air)	


¹⁾ Antireflection coating for visual observation and laser treatment:
450 - 810 nm

• Panfundus
Kontaktglas RETINA 145 L
Verre de contact RETINA 145 L
Contact glass RETINA 145 L

Optische Konfiguration des Kontaktglases RETINA 145 L
Configuration optique du verre de contact RETINA 145 L
Optical configuration of the contact glass RETINA 145 L

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4 Operating the product

The prerequisites for a successful application of contact glasses are an efficient anesthesia of the cornea and conjunctiva, which is not damaging the cornea epithelium.

The space between the eyes and the contact glass is to be filled with drops of 2 % Methocel, being an isotonic fluid non-damaging to the cornea and conjunctiva.

Fluctuations in temperature only exercise a minimal influence over the viscosity of the contact fluid. The fluid has the property of not seeping away when applied and the appearance of air bubbles is also extensively prevented.

Applying the contact glass

The contact glass can be applied to the anesthetized eye without any problem whilst the head of the patient remains supported on the headrest. The patient is requested to look upwards and the examiner draws the lower lid down and away from the eye.

If required, the upper lid can be drawn upwards somewhat with the thumb of the other hand. The contact glass is moistened with a drop of Methocel and placed on the conjunctiva and the patient then requested to look straight to the front, whereby the contact glass is permitted to glide over the cornea. Any air bubbles can be dispersed by tilting and turning the contact glass.

CAUTION

It is imperative to read the chapter "Operating safety" before using the product and to observe its precautions.

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Appendix A Maintenance

A.1 Cleaning and disinfecting of contact glasses

• Contact glasses are not dispatched in a disinfected state and should always be disinfected before use.

• Thorough cleaning under running cold tapwater during 30 to 60 seconds is mandatory before undertaking disinfecting routines. For more intensive cleaning in case of strong contamination carefully rub the contact glass with a cotton-wool ball under running water. In addition only mild soap or ether may be used.

• For disinfecting purposes we recommend an aqueous "Pantasept®" solution of 3 % during 10 minutes. Refer to the latest versions of the manufacturer's instructions and sales literature for details of the exact concentration, effecting time and application of disinfectant.

Alternative agents: only those products which are listed on HAAG-STREIT AG internet homepage (www.haag-streit.com) may be used.

• After disinfection the contact glasses must be rinsed thoroughly at least 10 to 60 minutes under running cold drinking water. Dry them with a clean, soft one-way tissue and store them in a clean, dry container.

• Only contact glasses, free of any residue, should be used. Residue may cause provoking or caustic reactions on the patient's eye.

A.2 Accessories for disinfecting and storing

For efficient disinfecting routines we recommend our "DESINSET", which comprises

- cleaning- and
- storage containers
- "Pantasept®" disinfecting solution

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Not allowed! Contact glasses can get damaged because of:

- Alcohol
- Acetone
- UV radiation
- Sterilization
- Immersing in fluid for more than 1 hour
- Temperature above 60 °C

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
(see also: [Appendix A](#) and [www.haag-streit.com](#))

30 - 60 seconds

Thorough cleaning under **running (cold) water** - removes Methocel and lacrimal fluid. In addition only mild soap and ether are allowed.

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reinigen / nettoyer / clean

Disinfect contact glasses:


see instructions of specific disinfectant. With "Pantasept®" HAAG-STREIT recommends:

10 minutes

in aqueous solution of 3 %

- No alcohol, no acetone!
- UV disinfecting and sterilization are not allowed!
- max. 60° C

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
min. 10 - max. 60 minutes

Clean thoroughly under **running (cold) drinking water**

3


spülen / rincer / rinse

Dry with a clean, **soft one-way tissue**


4

trocknen / sécher / dry

Place into **clean, dry storage container**

5

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Recommendation Use the HAAG-STREIT "DESINSET" (Accessories for Disinfection)

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Appendix B Technical specifications

Modell Modèle Model	Mit skleralem Rand Avec bord scléral With scleral flange	Laser	Spiegel Miroir Mirror	Bereich Aire Area	Vergrösserung Grossissement Magnification	Sehfeld Champ de vision Vision field	Seite / page / page	Spiegel Winkel Angle de miroir Mirror angle	r Kugel [mm] r Boule [mm] r Ball [mm]	ø Kontakt [mm] ø Contacte [mm] ø Contact [mm]	ø Kontakt aussen [mm] ø Contact externe [mm] ø Contact exterior [mm]	Höhe [mm] Hauteur [mm] Height [mm]	Gewicht [g] Poids [g] Weight [g]	Argon / Diode / YAG 532 / YAG 450 nm - 650 nm / 1064 nm	Argon / Diode / YAG 532 450 nm - 650 nm	Argon / Diode / YAG 532 450 nm - 810 nm	Stery Cup
630	.	3	Fundus / iridocornea (short)	0.91 x	.	22/23		59°/66°/73°	7.4	12	18	25.3	10.6	.	.	.	.
630 L	.	✓ 3	Fundus / iridocornea (short)	0.91 x	.	30/31		59°/66°/73°	7.4	12	18	25.8	11	✓	.	.	.
901	.	direct	Fundus	0.96 x	30°	12/13		.	7.4	12	15.5	12.8	1.5	.	.	.	.
902	.	1	Gonio	.	.	14/15		62°	7.4	12	15.5	21.7	4.2	.	.	.	.
902 S	✓	1	Gonio	.	.	14/15		62°	7.4	12	20	24	5.8	.	.	.	.
903	.	3	Fundus / iridocornea	0.91 x	.	20/21		59°/66°/73°	7.4	12	18.3	32	19	.	.	.	✓
903 L	.	✓ 3	Fundus / iridocornea	0.91 x	.	30/31		59°/66°/73°	7.4	12	18.3	32.5	20.2	✓	.	.	✓
903 S	✓	3	Fundus / iridocornea	0.91 x	.	20/21		59°/66°/73°	7.4	12	20.5	33.5	19.6	.	.	.	.
904	.	1	Ora serrata	.	.	18/19		62°	7.4	12	22	32.6	35.2	.	.	.	.
905	.	2	Gonio	.	.	16/17		62°/62°	7.4	12	15.5	21.7	5	.	.	.	.
905 S	✓	2	Gonio	.	.	16/17		62°/62°	7.4	12	20	24	5.8	.	.	.	.
906	.	3	Fundus / iridocornea (Baby)	.	.	22/23		59°/66°/73°	7	10	15.4	28	12.2	.	.	.	.
906 L	.	✓ 3	Fundus / iridocornea (Baby)	.	.	30/31		59°/66°/73°	7	10	15.4	28.5	12.8	✓	.	.	.
906 S	✓	3	Fundus / iridocornea (Baby)	.	.	22/23		59°/66°/73°	7	10	17.6	29.5	12.6	.	.	.	.
907	.	3	Fundus / iridocornea (Children)	.	.	22/23		59°/66°/73°	7.3	11	16.8	30.2	15.8	.	.	.	.
907 L	.	✓ 3	Fundus / iridocornea (Children)	.	.	30/31		59°/66°/73°	7.3	11	16.8	30.8	16.6	✓	.	.	.
1210	.	direct	Endothelium	2.2 x	.	24-27		.	7.3	12	18	25.5	7.8	.	.	.	.
CGAL	.	✓ 1	Gonio	1.5 x	.	28/29		58°	7.4	12	16	24	17.2	✓	.	.	.
CGIL	.	✓ direct	Iris	1.6 x	.	34/35		.	7.4	12	14	24	14.4	✓	.	.	.
CGPL	.	✓ direct	Pupillary region	1.5 x	.	38/39		.	7.4	12	15	13	4.4	✓	.	.	.
CGRL	.	✓ direct	Retina	0.75 x	64°	32/33		.	7.4	12	15.5	13	5.4	.	✓	.	.
CGVL	.	✓ direct	Vitreous	1.4 x	.	36/37		.	7.4	12	16	13	6.2	✓	.	.	.
RETINA 145 L	.	✓ direct	Panfundus	- 0.7 x	145°	40/41		.	7.4	12	18	27.5	26	.	.	✓	✓

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CE **HAAG-STREIT AG**